

R J Hill Laboratories Limited 1 Clyde Street Private Bag 3205 Hamilton 3240, New Zealand

Tel +64 7 858 2000 Fax +64 7 858 2001 Email mail@hill-labs.co.nz Web www.hill-labs.co.nz

ANALYSIS REPORT

Page 1 of 2

SPv1

Client: Greater Wellington
Contact: Dr John Drewry

C/- Greater Wellington
The Regional Council

PO Box 41 Masterton 5840 **Lab No:** 1694235 **Date Received:** 09-Dec-2016

Date Reported: 14-Dec-2016

Quote No: 82561 **Order No:** 228037

Client Reference:

Submitted By: Dr John Drewry

Sample Type: Soil								
	Sample Name:	Composite of 1, 2, 3 and 4						
	Lab Number:	1694235.5						
Total Recoverable Arsenic	mg/kg dry wt	22	-	-	-	-		
Total Recoverable Chromium	mg/kg dry wt	20	-	-	-	-		
Total Recoverable Copper	mg/kg dry wt	55	-	-	-	-		
Total Recoverable Iron	mg/kg dry wt	17,400	-	-	-	-		
Total Recoverable Lead	mg/kg dry wt	121	-	-	-	-		
Total Recoverable Manganese	e mg/kg dry wt	830	-	-	-	-		
Total Recoverable Zinc	mg/kg dry wt	300	-	-	-	-		
pH*	pH Units	7.1	-	-	-	-		

SUMMARY OF METHODS

The following table(s) gives a brief description of the methods used to conduct the analyses for this job. The detection limits given below are those attainable in a relatively clean matrix. Detection limits may be higher for individual samples should insufficient sample be available, or if the matrix requires that dilutions be performed during analysis.

Sample Type: Soil							
Test	Method Description	Default Detection Limit	Sample No				
Environmental Solids Sample Preparation	Air dried at 35°C and sieved, <2mm fraction. Used for sample preparation. May contain a residual moisture content of 2-5%.	-	1-5				
Soil Prep Dry & Sieve for Agriculture	Air dried at 35°C and sieved, <2mm fraction.	-	5				
Total Recoverable digestion	Nitric / hydrochloric acid digestion. US EPA 200.2.	-	5				
Composite Environmental Solid Samples*	Individual sample fractions mixed together to form a composite fraction.	-	1-4				
Total Recoverable Arsenic	Dried sample, sieved as specified (if required). Nitric/Hydrochloric acid digestion, ICP-MS, screen level. US EPA 200.2.	2 mg/kg dry wt	5				
Total Recoverable Chromium	Dried sample, sieved as specified (if required). Nitric/Hydrochloric acid digestion, ICP-MS, screen level. US EPA 200.2.	2 mg/kg dry wt	5				
Total Recoverable Copper	Dried sample, sieved as specified (if required). Nitric/Hydrochloric acid digestion, ICP-MS, screen level. US EPA 200.2.	2 mg/kg dry wt	5				
Total Recoverable Iron	Dried sample, sieved as specified (if required). Nitric/Hydrochloric acid digestion, ICP-MS, screen level. US EPA 200.2.	40 mg/kg dry wt	5				
Total Recoverable Lead	Dried sample, sieved as specified (if required). Nitric/Hydrochloric acid digestion, ICP-MS, screen level. US EPA 200.2.	0.4 mg/kg dry wt	5				
Total Recoverable Manganese	Dried sample, sieved as specified (if required). Nitric/Hydrochloric acid digestion, ICP-MS, screen level. US EPA 200.2.	1.0 mg/kg dry wt	5				
Total Recoverable Zinc	Dried sample, sieved as specified (if required). Nitric/Hydrochloric acid digestion, ICP-MS, screen level. US EPA 200.2.	4 mg/kg dry wt	5				
рН*	1:2 (v/v) soil : water slurry followed by potentiometric determination of pH.	0.1 pH Units	5				



These samples were collected by yourselves (or your agent) and analysed as received at the laboratory.

Samples are held at the laboratory after reporting for a length of time depending on the preservation used and the stability of the analytes being tested. Once the storage period is completed the samples are discarded unless otherwise advised by the client.

This report must not be reproduced, except in full, without the written consent of the signatory.

Ara Heron BSc (Tech)

Client Services Manager - Environmental